

Product Reference Guide

Polycarbonate

Photochromic Gray and Brown

Aspheric Finished Single Vision

64mm Aspheric Finished Single Vision

Sphere Power Up to -2.00 cyl	Base Curve	Curve Radius	Center Thickness
+0.25	4.28	123.832	1.6
+0.50	4.28	123.832	1.7
+0.75	4.28	123.832	1.8
+1.00	4.73	112.051	2.0
+1.25	4.73	112.051	2.2
+1.50	4.73	112.051	2.5
+1.75	5.16	102.713	2.6
+2.00	5.16	102.713	2.8

Base Curve Selection

	Cylinders								
Spheres	0.00	-0.25	-0.50	-0.75	-1.00	-1.25	-1.50	-1.75	-2.00
2.00									
1.00					64mm				
.25									

Product Reference Guide

Polycarbonate

Photochromic Gray and Brown

Aspheric Finished Single Vision

75mm Aspheric Finished Single Vision

Sphere Power Up to -1.50 cyl	Base Curve	Curve Radius	Center Thickness
Plano	4.28	123.832	1.5
-0.25	4.28	123.832	1.5
-0.50	4.28	123.832	1.5
-0.75	4.28	123.832	1.5
-1.00	4.28	123.832	1.5
-1.25	3.83	138.381	1.5
-1.50	3.83	138.381	1.5
-1.75	3.83	138.381	1.5
-2.00	3.83	138.381	1.5
-2.25	3.38	156.805	1.5
-2.50	3.38	156.805	1.5
-2.75	3.38	156.805	1.5
-3.00	3.38	156.805	1.5
-3.25	2.93	180.887	1.4
-3.50	2.93	180.887	1.4
-3.75	2.93	180.887	1.4
-4.00	2.93	180.887	1.4
-4.25	2.48	213.710	1.4
-4.50	2.48	213.710	1.4
-4.75	2.48	213.710	1.3
-5.00	2.48	213.710	1.3
-5.25	1.81	292.818	1.3
-5.50	1.81	292.818	1.3
-5.75	1.81	292.818	1.3
-6.00	1.81	292.818	1.3

