

WHY CLEAR BLUE FILTER™?

Clear Protection
from Damaging
Light Rays



All eyes are vulnerable.

From sunlight to digital devices, ultraviolet and high-energy blue light is everywhere. Clear Blue Filter™ lenses may help protect eyes from the harmful effects of powerful rays.

WHY?

Help safeguard eyes from cataracts, macular degeneration and other serious conditions linked to UV and high-energy blue light.^{1,2} More exposure means more risk, especially for kids.



250% increase
in LED bulb
sales³

70% more high-energy
blue light transmits
to a child's retina⁴

¹ "The Family of Sunlight-Related Eye Diseases" Optometry & Vision Science, Volume 71 | ² "Is sunlight exposure a risk factor for age-related macular degeneration?" British Journal of Ophthalmology, Volume 97 | ³ "GE Lighting plans to stop CFL bulb production by years end." Feb 2016. <http://www.fastcodesign.com> | ⁴ Sidney Lerman, M.D. "Clinical Light Damage to the Eye"

Provide the best protection.

Block UVB/UVA rays and filter the most damaging, highest-energy blue light with Clear Blue Filter™

HOW?

With Clear Blue Filter™, safety is built into the lens, not the coating.

Designed for daily wear.

These lenses are thin, lightweight, impact-resistant and affordable—making them perfect for everyday use.



**SKIN CANCER FOUNDATION
SEAL OF APPROVAL**

Recommended as an effective UV filter
for the eyes and surrounding skin.





It's the clear choice.

Give them crystal clear lenses, more comfort from screen time and eye protection all in one with Clear Blue Filter™.

96% of consumers rate lens clarity as highly important

3x more blue light filtered than standard polycarbonate at 420nm

100% UV blocking up to 400nm

Availability

| Style | Material | Diameter (mm) | Power Range | Bases | CYL/Add |
|----------------------------|-----------------|---------------------------|-----------------|--|--------------------|
| SFSV HC | Poly | 75mm | | 0.50, 1.00, 2.00, 3.00, 4.00, 5.25, 6.25, 7.50, 8.50, 9.75 | |
| FSV HC | Poly | Spherical 65mm, 75mm | +4.00 to -8.00 | | out to a -2.00 cyl |
| A FSV Super-Hydro Vivid AR | Poly | Aspheric 65mm, 75mm | +4.00 to -6.00 | | out to a -2.00 cyl |
| SFSV HC | Plastic 1.50 | 75mm | | 0.50, 1.25, 2.25, 3.25, 4.25, 5.25, 6.25, 7.25, 8.25, 9.25 | |
| FSV HC | Plastic 1.50 | 70mm | +4.00 to -6.00 | | out to a -2.00 cyl |
| SFSV HC | High Index 1.67 | 60mm, 70mm, 75mm | | 0.50, 1.00, 2.50, 4.00, 5.00, 6.00, 7.00, 8.00, 9.00, 10.00, 11.00 | |
| A FSV Super-Hydro Vivid AR | High Index 1.67 | Aspheric 65mm, 70mm, 75mm | +6.00 to -12.00 | | out to a -2.00 cyl |

READY TO ORDER?

Rx Orders
Phone:1-800-328-3449
VelRX@hoya.com

Add Clear Blue Filter™ lenses to your portfolio today.
[visionease.com](https://www.visionease.com)

Stock Orders
Phone: 1-833-317-0410
orders@visionease.com