

Vision Ease Stock Lens Availability Guide 2025

CLEAR PLASTIC 1.50						
Material	Product	Diameter	Power Range	Colors	Bases	Cyl/Adds
FINISHED SINGLE VISION						
Plastic 1.50	VISION EASE FSV Uncoated	70mm	+4.00 to -6.00			out to a -2.00 cyl
	VISION EASE FSV Hardcoat	70mm	+4.00 to -6.00			out to a -2.00 cyl
	VISION EASE FSV AR	70mm	+4.00 to -6.00			out to a -2.00 cyl
SEMI-FINISHED SINGLE VISION						
Plastic 1.50	VISION EASE SFSV Hardcoat	70mm			Plano, 2.25, 4.25, 5.25, 6.25, 7.25, 8.25, 10.25	
	VISION EASE SFSV Hardcoat	75mm			1.25, 2.25, 3.25, 4.25, 5.25, 6.25, 7.25, 8.25, 10.25	
	VISION EASE SFSV Uncoated	65mm			Plano, 2.25, 4.25, 5.25, 6.25, 8.25, 10.25	
	VISION EASE SFSV Uncoated	70mm			Plano, 2.25, 4.25, 5.25, 6.25, 7.25, 8.25, 10.25	
	VISION EASE SFSV Uncoated Extra Thick	70mm			Plano, 2.25, 4.25, 6.25, 8.25, 10.25	
	VISION EASE SFSV Uncoated	75mm			1.25, 2.25, 3.25, 4.25, 5.25, 6.25, 7.25, 8.25, 10.25	
	VISION EASE SFSV Uncoated	80mm			2.25, 4.25, 5.25, 6.25, 7.25, 8.25	
SEMI-FINISHED BIFOCALS AND TRIFOCALS						
Plastic 1.50	VISION EASE D28 Hardcoat	70mm			0.50, 2.25, 4.25, 6.25, 8.25, 10.25	+1.00 to +3.00
	VISION EASE D35 Hardcoat	70mm			0.50, 2.25, 4.25, 6.25, 8.25, 10.25	+1.00 to +3.00
	VISION EASE RD24 Hardcoat	70mm			4.25, 6.25, 8.25	+1.00 to +3.00
	VISION EASE 7x28 Hardcoat	70mm			2.25, 4.25, 6.25, 7.25, 8.25	+1.50 to +3.00
	VISION EASE D25 Uncoated	70mm			0.50, 2.25, 4.25, 6.25, 8.25, 10.25	+0.75 to +3.50
	VISION EASE D28 Uncoated	70mm			0.50, 2.25, 4.25, 6.25, 7.25, 8.25, 10.25	+0.75 to +4.00
	VISION EASE D28 Uncoated	75mm			4.25, 6.25, 8.25	+0.50 to +4.00
	VISION EASE D28 Uncoated	80mm			4.25, 6.25, 8.25	+1.00 to +3.00
	VISION EASE D35 Uncoated	70mm			2.25, 4.25, 6.25, 8.25, 10.25	+0.75 to +4.00
	VISION EASE D35 Uncoated	75mm			4.25, 6.25, 8.25	+1.00 to +4.00
	VISION EASE RD24 Uncoated	70mm			0.50, 2.25, 4.25, 6.25, 8.25, 10.25	+1.00 to +3.00
	VISION EASE 7x28 Uncoated	70mm			2.25, 4.25, 6.25, 7.25, 8.25	+1.50 to +3.00
	VISION EASE 7x28 Uncoated	75mm			4.25, 6.25, 8.25	+1.50 to +3.00
	VISION EASE 8x35 Uncoated	70mm			4.25, 6.25, 8.25	+1.50 to +3.00
	VISION EASE 8x35 Uncoated	75mm			4.25, 6.25, 8.25	+1.50 to +3.00
CLEAR HIGH INDEX 1.67						
Material	Product	Diameter	Power Range	Colors	Bases	Cyl/Adds
FINISHED SINGLE VISION						
High Index 1.67	VISION EASE Aspheric FSV SHAR	65mm, 70mm, 75mm	+6.00 to -12.00			out to a -2.00 cyl
SEMI-FINISHED SINGLE VISION						
High Index 1.67	VISION EASE Aspheric SFSV Hardcoat	70mm, 75mm			0.75, 2.00, 3.25, 5.00, 6.00, 7.00, 8.00, 10.00	
	VISION EASE SFSV (Coated/Uncoated, MR7/MR10)	70mm, 75mm			0.75, 2.00, 3.50, 4.75, 6.00, 7.50, 9.00	

CLEAR POLYCARBONATE						
Material	Product	Diameter	Power Range	Colors	Bases	Cyl/Adds
FINISHED SINGLE VISION						
Polycarbonate	VISION EASE FSV Advanta Hardcoat	65mm, 75mm	+4.00 to -8.00			out to a -2.00 cyl
	VISION EASE FSV Hardcoat	65mm, 73mm	+4.00 to -8.00			out to a -2.00 cyl
	VISION EASE FSV Hardcoat (High Cylinder)	65mm, 73mm	+4.00 to -8.00			out to a -2.00 cyl
	VISION EASE Aspheric FSV Tintable Hardcoat	65mm	+4.00 to +0.25			out to a -2.00 cyl
	VISION EASE FSV Tintable Hardcoat	75mm	Plano to -8.00			out to a -2.00 cyl
	VISION EASE FSV Safety Continua Hardcoat	70mm, 75mm	+4.00 to -8.00			out to a -2.00 cyl
	VISION EASE FSV VELA AR	65mm, 73mm	+4.00 to -6.00			out to a -4.00 cyl
	VISION EASE FSV VELA AR (High Cylinder)	Spherical 65mm plus, 73mm minus	+4.00 to -6.00			out to a -4.00 cyl
	VISION EASE FSV VELA Super-Hydro AR	Spherical 65mm plus, 73mm minus	+4.00 to -6.00			out to a -4.00 cyl
	VISION EASE FSV VELA Super-Hydro AR (High Cylinder)	Spherical 65mm plus, 73mm minus	+4.00 to -6.00			out to a -4.00 cyl
	VISION EASE Aspheric FSV Tegra Hardcoat	65mm, 75mm	+4.00 to -6.00			out to a -2.00 cyl
	VISION EASE Aspheric FSV Premium SHAR	Aspheric 64mm plus, 75mm minus	+4.00 to -6.00			out to a -2.00 cyl
SEMI-FINISHED SINGLE VISION						
Polycarbonate	VISION EASE SFSV Hardcoat	71mm, 75mm, 80mm			0.50, 1.00, 2.00, 3.00, 4.00, 5.25, 6.25, 7.50, 8.50, 9.75	
	VISION EASE SFSV Tegra	75mm			0.50, 1.00, 2.00, 3.00, 4.00, 5.25, 6.25, 7.50, 8.50, 9.75	
	VISION EASE Aspheric SFSV Tegra	75mm			0.50, 1.00, 2.00, 3.00, 4.00, 5.25, 6.25, 7.50, 8.50	
	VISION EASE SFSV Strippable Hardcoat	75mm			0.50, 1.00, 2.00, 3.00, 4.00, 5.25, 6.25, 7.50, 8.50, 9.75	
SEMI-FINISHED PROGRESSIVE						
Polycarbonate	VISION EASE SF Novel Decentered Progressive	80mm			1.50, 3.00, 4.00, 5.00, 7.00	+0.75 to +4.00
SEMI-FINISHED BIFOCALS & TRIFOCALS						
Polycarbonate	VISION EASE D28 Bifocal	75mm	+6.00 to -11.00		2.00, 4.00, 6.25, 8.50	+1.00 to +3.00
	VISION EASE Tegra D28 Bifocal	75mm	+6.00 to -9.00		2.50, 3.25, 4.25, 5.50, 6.50	+1.00 to +3.00
	VISION EASE D35 Bifocal	75mm	+6.00 to -9.00		2.00, 4.00, 6.25, 8.50	+1.00 to +3.00
	VISION EASE 7x28 Trifocal	75mm	+6.00 to -9.00		2.00, 4.00, 6.25, 8.50	+1.50 to +4.00
	VISION EASE 8x35 Trifocal	75mm	+6.00 to -9.00		2.00, 4.00, 6.25, 8.50	+1.50 to +3.00

BLUE LIGHT PROTECTION						
Material	Product	Diameter	Power Range	Colors	Bases	Cyl/Adds
FINISHED SINGLE VISION						
Plastic 1.50	VISION EASE Clear Blue Filter FSV Hardcoat	70mm	+4.00 to -6.00			out to a -2.00 cyl
	VISION EASE Clear Blue Filter FSV AR	70mm	+4.00 to -6.00			out to a -2.00 cyl
Polycarbonate	VISION EASE Clear Blue Filter Aspheric FSV Super-Hydro Vivid AR	Aspheric 65mm, 75mm	+4.00 to -6.00			out to a -2.00 cyl
	VISION EASE Clear Blue Filter FSV Spherical Hardcoat	Spherical 65mm, 75mm	+4.00 to -8.00			out to a -2.00 cyl
High Index 1.67	VISION EASE Clear Blue Filter Aspheric FSV Super-Hydro Vivid AR	Aspheric 65mm, 70mm, 75mm	+6.00 to -12.00			out to a -2.00cyl
SEMI-FINISHED SINGLE VISION						
Plastic 1.50	VISION EASE Clear Blue Filter SFSV Hardcoat	75mm			0.50, 1.25, 2.25, 3.25, 4.25, 5.25, 6.25, 7.25, 8.25, 9.25	
	VISION EASE Clear Blue Filter SFSV Uncoated	75mm			0.50, 1.25, 2.25, 3.25, 4.25, 5.25, 6.25, 7.25, 8.25, 9.25	
Polycarbonate	VISION EASE Clear Blue Filter SFSV Hardcoat	75mm			0.50, 1.00, 2.00, 3.00, 4.00, 5.25, 6.25, 7.50, 8.50, 9.75	
High Index 1.67	VISION EASE Clear Blue Filter SFSV Hardcoat	65mm, 70mm, 75mm			0.50, 1.00, 2.50, 4.00, 5.00, 6.00, 7.00, 8.00, 9.00, 10.00, 11.00	
SEMI-FINISHED BIFOCALS AND TRIFOCALS						
Plastic 1.50	VISION EASE Clear Blue Filter D28 Uncoated	70mm			0.50, 2.00, 4.00, 6.00, 8.00	+1.50 to +3.50
	VISION EASE Clear Blue Filter D35 Uncoated	70mm			2.00, 4.00, 6.00, 8.00	+1.00 to +3.50
	VISION EASE Clear Blue Filter 7x28 Uncoated	70mm			4.00, 6.00, 8.00	+1.50 to +3.50
	VISION EASE Clear Blue Filter 8x35 Uncoated	70mm			2.00, 4.00, 6.00, 8.00	+1.50 to +3.00
Polycarbonate	VISION EASE Clear Blue Filter D28 Hardcoat	75mm			0.50, 2.00, 4.00, 6.25, 8.50	+0.75 to +4.00
	VISION EASE Clear Blue Filter D35 Hardcoat	75mm			2.00, 4.00, 6.25, 8.50	+1.00 to +3.00
	VISION EASE Clear Blue Filter 7x28 Hardcoat	75mm			2.00, 4.00, 6.25, 8.50	+1.50 to +4.00
High Index 1.67	VISION EASE Clear Blue Filter D28 Uncoated	75mm			1.00, 3.00, 5.00, 7.00	+1.00 to +3.50

1 Double Hardcoat (direct apply AR/strip)

2 Single Hardcoat (optional over-dip)

*Brown power range = +2.00 to -4.00

POLARIZED						
Material	Product	Diameter	Power Range	Colors	Bases	Cyl/Adds
FINISHED SINGLE VISION						
Plastic 1.50	VISION EASE Polarized Plano FSV Coated	70mm	Plano	Gray 3, Brown 3		out to a -2.00 cyl
Polycarbonate	VISION EASE Polarized Plano FSV	75mm	Plano; Base 6.00, 8.00	Gray, Brown		out to a -2.00 cyl
SEMI-FINISHED SINGLE VISION						
Plastic 1.50	VISION EASE Polarized Plano SFSV Uncoated	76mm		Gray 3, Brown 3	6.25	
	VISION EASE Polarized SFSV Coated and Uncoated	77mm		Gray 3, Brown 3	2.25, 3.25, 4.25, 5.25, 6.25, 7.25, 8.25	
Phoenix	HOYA Polarized SFSV Coated	74mm		Gray, Brown	2.25, 4.25, 6.25, 8.25	
Polycarbonate	COPPERTONE SFSV	75mm		Gray, Brown, Green	0.50, 1.00, 2.00, 3.00, 4.00, 5.25, 6.25, 7.50, 8.50	
	COPPERTONE SFSV	80mm		Gray, Brown, Green	4.00, 6.25, 8.50	
	COPPERTONE Plano SFSV	75mm		Gray, Brown, Green	6.00, 8.00	
	VISION EASE Polarized SFSV	75mm		Gray 3, Brown 3, Classic Green*, Blue**	0.50, 1.00, 2.00, 3.00, 4.00, 5.25, 6.25, 7.50, 8.50	
	VISION EASE Polarized SFSV	80mm		Gray 3	4.00, 6.25, 8.50	
	VISION EASE Polarized Mirror Coat SFSV	75mm		Silver, Gold, Blue (Polar Gray base)	2.00, 4.00, 6.25, 7.50	
	VISION EASE Polarized Plano SFSV	75mm		Gray 3, Brown 3	6.00, 8.00	
	VISION EASE Polarized Tints SFSV	75mm		Gray 2	2.00, 4.00, 6.25, 8.50	
SEMI-FINISHED PROGRESSIVES						
Polycarbonate	COPPERTONE Novel Decentered Progressive Hardcoat	80mm, 85mm		Gray, Brown, Green	1.50, 3.00, 4.00, 5.00, 7.00	+1.00 to +3.00
	COPPERTONE Novel Decentered Progressive Hardcoat	75mm	Plano	Gray, Brown		
	VISION EASE Polarized Novel Decentered Progressive Hardcoat	80mm, 85mm		Gray 3, Brown 3	1.50, 3.00, 4.00, 5.00, 7.00	+1.00 to +3.00
SEMI-FINISHED BIFOCALS AND TRIFOCALS						
Plastic 1.50	VISION EASE Polarized D28 Bifocal Coated and Uncoated	77mm		Gray 3, Brown 3	2.25, 4.25, 6.25, 8.25	+1.00 to +3.00
Polycarbonate	COPPERTONE D28 Bifocal	75mm		Gray, Brown, Green	2.00, 4.00, 6.25, 8.50	+1.00 to +3.00
	COPPERTONE 7x28 Trifocal	75mm		Gray, Brown	4.00, 6.25	+1.00 to +3.00
	VISION EASE Polarized D28 Bifocal	75mm		Gray 3, Brown 3, Classic Green***	2.00, 4.00, 6.25	+1.00 to +3.00
	VISION EASE Polarized 7x28 Trifocal	75mm		Gray 3, Brown 3, Classic Green	2.00, 4.00, 6.25	+1.50 to +3.00

*No 8.50 Base, **4.00, 6.25, 8.50 Base, ***+1.25 to +3.00 add, No 2.00 Base



PHOTOCHROMIC						
Material	Product	Diameter	Power Range	Colors	Bases	Cyl/Adds
FINISHED SINGLE VISION						
Polycarbonate	VISION EASE Photochromic Aspheric FSV	65mm, 75mm	+2.00 to -6.00*	Gray, Brown		out to a -2.00 cyl
	VISION EASE Photochromic Aspheric FSV SHAR	65mm, 75mm	+2.00 to -6.00*	Gray, Brown		out to a -2.00 cyl
SEMI-FINISHED SINGLE VISION						
Polycarbonate	VISION EASE Photochromic SFSV	75mm		Gray, Brown	0.50, 1.00, 2.00, 3.00, 4.00, 5.25, 6.25, 7.50	
SEMI-FINISHED PROGRESSIVE						
Polycarbonate	VISION EASE Photochromic Novel Decentered Progressive	80mm		Gray	1.50, 3.00, 4.00, 5.00, 7.00	+1.00 to +3.00
SEMI-FINISHED BIFOCALS AND TRIFOCALS						
Plastic 1.50	LENSCHOICE Photo D28 Hardcoat	75mm		Gray, Brown	2.00, 4.00, 6.00, 8.00	+1.00 to +3.00
	LENSCHOICE Photo D35 Hardcoat/Uncoated	75mm		Gray, Brown	2.25, 4.25, 6.25, 8.25	+1.00 to +3.50
	LENSCHOICE Photo 7x28 Trifocal Coated	70mm, 75mm		Gray, Brown	2.25, 4.25, 6.25, 8.25	+1.50 to +3.00
	LENSCHOICE Photo 8x35 Trifocal Coated	75mm		Gray, Brown	4.25, 6.25, 8.25	+1.50 to +3.00
Polycarbonate	VISION EASE Photochromic D28 Bifocal	75mm		Gray, Brown	0.50, 2.00, 4.00, 6.25, 8.50	+1.00 to +3.00
	VISION EASE Photochromic D35 Bifocal	75mm		Gray, Brown	2.00, 4.00, 6.25	+1.00 to +3.00
	VISION EASE Photochromic 7x28 Trifocal	75mm		Gray, Brown	2.00, 4.00, 6.25	+1.50 to +3.00